- GB

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
 **EC REP**
Donawa Lifescience Consulting Srl
Piazza Albania, 10
00153 Rome, Italy

GB

ONGUARD® Closed System Transfer Device (CSTD) containing TEVADAPTOR® Drug Reconstitution and Administration Components

INDICATIONS FOR USE:
TEVADAPTOR is a Closed System Drug Transfer Device (CSTD) that mechanically prohibits the release of the drug in vapor, aerosol or liquid form during preparation and administration, and prevents the introduction of microbial and airborne contaminants into the drug or fluid path, allowing the system to minimize exposure of individuals, healthcare personnel, and the environment to hazardous drugs.

- NOTES:**
- Use Aseptic technique.
 - To expel air from syringe, draw half of required volume, tap inverted syringe lightly, gently depress plunger to expel air only, then draw on plunger to total required volume.
 - Replace tubing sets per CDC guidelines and/or institutional protocol.
 - Disposal of components should be according to local regulations and procedures for handling of hazardous drugs.
 - To maintain closed system, do not disconnect tubing sets.
 - The various component parts of the TEVADAPTOR system are available separately if required.
 - Refer to B. Braun Medical Inc. for additional training information.
 - TEVADAPTOR has been shown to prevent microbial ingress when the septa of the TEVADAPTOR components are pierced up to 10 times over a 72 hour period
 - This device is not intended to be a substitute for the drug reconstitution process when sterile compounding is required, which requires ISO class 5 environment admixing or the use of a Biological Safety Cabinet (BSC) or laminar flow hood.

- PRECAUTIONS:**
- Do not fill vial in an inverted position.
 - Do not disconnect *Vial Adaptor/Syringe Adaptor* from *Vial/Syringe* once connected.
 - Do not use with 13 mm vials of lyophilized drugs.
 - The *Syringe Adaptor/Syringe Adaptor Lock* hub is made of polycarbonate which is known to be incompatible with certain drugs. Do not use TEVADAPTOR products with drugs containing the solvent N,N-dimethylacetamide (e.g. Treanda, busulfan and amsacrine). Refer to drug manufacturer's package insert regarding drug preparation and compatibility.
 - The fluid path of the TEVADAPTOR system contains the following materials: ABS, Polyisoprene, Polycarbonate, and DEHP-free PVC. Refer to drug manufacturer's package insert regarding drug preparation and compatibility.
 - To prevent microbial ingress, swab the septa of TEVADAPTOR® components vigorously with 70% IPA for 5 seconds before each use.
 - When using the TEVADAPTOR system the user should follow the standard precautions which should be taken when handling hazardous drugs with drug transfer systems.
 - *Vial Adaptor* can be used up to 72 hours, however, storage of medication with the TEVADAPTOR system should not exceed the duration defined in the medication label.
 - Closed gravity administration using the TEVADAPTOR syringe adaptor allows a maximum saline flow rate of 90 mL/minute.

- Instructions for Use:**
1. Remove protective cap from vial **A**
 2. While keeping the vial on a flat surface, **A** insert the spike of the *Vial Adaptor* in the center of the rubber surface. **C**
 - For 13 mm vials, use the *Vial Converter* first by inverting the vial and pushing into converter. **B** *

- For 20 mm vials, use 20 mm *Vial Adaptor*
 - For 28 mm vials, use 28 mm *Vial Adaptor*
3. Connect *Syringe Adaptor/Syringe Adaptor Lock* to a luer lock syringe. **D**
 - Only legally marketed sterile, single use, syringes having female luer lock connectors meeting the requirements of ISO 594-2 are to be used with the *Syringe Adaptor/Syringe Adaptor Lock*.
 4. Select one of the TEVADAPTOR *Secondary Sets*. Close tubing clamp. Spike set into bag spike port. **E**
 5. To prime *Connecting Set*: (Do not remove tubing cap.) Position fluid bag above tubing set, open tubing clamp, and allow fluid to prime to tubing cap. Reclose clamp.

- For reconstitution (powder for injection) and Drug Transfer**
1. This device is not intended to be a substitute for the drug reconstitution process when sterile compounding is required, which requires ISO class 5 environment admixing or the use of a Biological Safety Cabinet (BSC) or laminar flow hood.
 2. Remove *Vial Adaptor* protective cap **F** and remove *Syringe Adaptor/Syringe Adaptor Lock* cap. **F1**
 3. Draw desired amount of diluent from a vial **G** or a bag. **I**
 - a. Click a new *Vial Adaptor* onto diluent vial, connect *Syringe Adaptor/Syringe Adaptor Lock*, and draw diluent from vial **G**, or
 - b. Connect *Syringe Adaptor/Syringe Adaptor Lock* to *Secondary Set* and draw desired amount of diluent from bag. **J**
 4. Disconnect *Syringe Adaptor/Syringe Adaptor Lock* slowly by pinching wings. **H**
 5. Connect *Syringe Adaptor/Syringe Adaptor Lock* to *Vial Adaptor* on drug vial, inject diluent, and dissolve drug. **K**
 6. Invert and withdraw desired volume from vial into syringe. **G**
 7. Disconnect *Syringe Adaptor/Syringe Adaptor Lock* slowly from *Vial Adaptor* by pinching wings. **H**
(Partial vial: cap *Vial Adaptor* and store vial in perpendicular position only for further use) **I**
NOTE: Before recapping/transporting the *Syringe Adaptor/Syringe Adaptor Lock* and syringe containing the drug, swab the cap and Tevadaptor septa with 70% IPA for 5 seconds, place it on the port and transport per institutional protocol.

- For Fluid Drug Transfer (drugs in solution)**
1. Remove *Vial Adaptor* protective cap. **F**
 2. Connect *Syringe Adaptor/Syringe Adaptor Lock* to *Vial Adaptor*, invert and withdraw desired drug volume. **G**
 3. Disconnect *Syringe Adaptor/Syringe Adaptor Lock* slowly from *Vial Adaptor* by pinching wings. **H**
(Partial vial: cap *Vial Adaptor* and store vial for further use.) **I**

- For IV Push/Bolus Administration**
1. Remove the cap from the *Syringe Adaptor/Syringe Adaptor Lock*.
 2. Firmly connect *Luer Lock Adaptor* **M** to appropriate primed female luer access of patient IV line.
 3. Connect *Syringe Adaptor/Syringe Adaptor Lock* and administer.
 4. Disconnect *Syringe Adaptor/Syringe Adaptor Lock* slowly by pinching wings **H** and dispose.
NOTE: To maintain closed system, do not disconnect *Luer Lock Adaptor* from patient IV line.

- For IV Drip Infusion**
1. Connect *Syringe Adaptor/Syringe Adaptor Lock* to *Secondary Sets* and inject drug into the bag. **L**
 2. Rinse spike and *Syringe Adaptor/Syringe Adaptor Lock* by withdrawing 2 to 3 mL of solution from bag and reinject.
 3. Disconnect *Syringe Adaptor/Syringe Adaptor Lock* from *Secondary Sets*. Transport per institutional protocol.
 4. Attach *Connecting Set* to appropriate female luer access of patient IV line or insert spike of primary IV set into *Spike Port Adaptor*. Prime set if needed.

- Use of Syringe Adaptor at a distal end of IV set**
1. Firmly connect *Syringe Adaptor* to distal end of the set.
 2. Connect *Luer Lock Adaptor* **M** to appropriate female luer access site of patient IV line.
 3. Connect *Syringe Adaptor* to *Luer Lock Adaptor* and administer.
 4. Disconnect *Syringe Adaptor* slowly by pinching wings and dispose of IV Set.
NOTE: To maintain closed system, do not disconnect *Luer Lock Adaptor* from patient IV line.

- TEVADAPTOR system components list:**
REF: 412111 - TEVADAPTOR 20 mm & 13 mm *Vial Adaptor* with TOXI-GUARD® Activated Charcoal Filter
REF: 412119 - TEVADAPTOR 28 mm *Vial Adaptor* with TOXI-GUARD Activated Charcoal Filter
For use in combination with *Syringe Adaptor*.
REF: 412118 - TEVADAPTOR Syringe Adaptor
For use in combination with *Vial Adaptor*, *Luer Lock Adaptor* and *Tevadaptor Sets*.
REF: 412126 - TEVADAPTOR Syringe Adaptor Lock
NOTE: Once the syringe is connected to Syringe Adaptor Lock, it cannot be disconnected.
REF: 412114 - TEVADAPTOR Luer Lock Adaptor
For use in combination with female luer lock connectors and *Syringe Adaptor*.

- TEVADAPTOR Secondary sets**
REF: 412112 - TEVADAPTOR Connecting Set
REF: 412115 - TEVADAPTOR Connecting Set with ULTRASITE® Injection Site
REF: 412113 - TEVADAPTOR Spike Port Adaptor
REF: 412116 - TEVADAPTOR Spike Port Adaptor with ULTRASITE Injection Site
REF: 412120 - TEVADAPTOR Secondary IV Administration Set with ULTRASITE Injection Site
For use in combination with *Syringe Adaptor*.

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ONGUARD is a registered trademark of B. Braun Medical Inc.
TEVADAPTOR and TOXI-GUARD are registered trademarks of Teva Medical Ltd.